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INSTRUCTIONS

When operating, insertion DC transformer into AC 220V socket, the another end should be inserted into the control box of health belt, at this time the display screen shows "OFF".

This microcomputer digital display control box has three button switches:

1. On/Off button: press this button to start this machine, and entering automatic circulating operation, completing one function every minute. After finishing one circulating massage within the first fixed 5 minutes, the machine will be off automatically, the display screen shows "OFF".

2. Function button: this microcomputer control box has 8 time programs. When starting, it will automatically set 5 minutes, later press this button, the fixed times displayed by the display screen are separately 5 minutes, 10 minutes, 15 minutes, 20 minutes, 25 minutes, 30 minutes, 35 minutes and 40 minutes, they can be selected freely.

After operating for a period of time, the user has sense of itch, paralysis or multi sweat, it is the normal appearance, it is caused by magical effect of over magnetic wave, vibration and massage.

Usually, it is suitable to operate for 15 minutes once, 2-4 time every day.

You should tightly bind the health belt on the part you want to thin, binding it tighter, the effect is better.

In order to improve metabolism better, always drink a cup of boiled water of fibre juice after operating.

When you are in a car, use power of igniter of the car to enjoy the above functions equally.

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POINTS FOR ATTENTION

This product uses 12V DC low voltage power, it is easy to operate, safety and reliable. If you don't use it, always cut off voltage stabilizing power and electric supply.

This belt is limited to use indoor and car, it should be avoided damp, and put in a dry place.

Don't use it after dinner within one hour. Pregnant woman, patients of heart disease, cardiovascular and cranial vascular disease are forbidden to use.

Don't disassemble and repair it by yourself if there is a breakdown. It should be repaired by Maintenance Department of the shop.

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BRIEFINTRODUCTION

Series of R.T-9803 is the new generation micro computer overmagnetic wave massage health belt, it is a new variety of the 21 st Century, adopting multi-powerful vibrating massage balls controlled by microcomputer, and displaying various massage form and time parameter with digital type, putting a permanent magnetic steel on relevant acupoint to produce powerful magnetic field. They can be wrapped on waist, neck, thigh and arm to make rolling, pinching, pushing, finger pressing and circulating massage. It can quicken blood circulation of human body, let the unnecessary fat to move directly, and taking the effect of separating to unnecessary fat gathered by intestinal mucosa of abdomen, quickly removing redundant fat, it will let fat and uncomfortable abdomen to become slim and fit.

It also can adjust organism organization, strengthen metabolism, clear and activate the channels and collaterals, remove fatigue, they will let you to keep vigorous energy, especially applying to ache diseases caused by damage of soft tissue of muscle of waist, fascia, ligament, joint bladder and etc., adjust circulation of qi(vital energy) and blood through transmission of main and collateral channels, and then reach calm and relieve pain, remove convulsion of muscle. After vibrating and massaging, it will, it really is the necessary, practical and best health goods for every household.

Connect With Us

Comfortable

Loosening

Reducing weight