

Multi-Lead Cables and Electrodes



Previous



Next

Table of Contents

Product Index

Search

Home


Contact Us



Previous



Next




5-Lead Clinical ECG Cables and Electrodes

A Set of 5 color-coded ECG Limb Cables, labelled RA, LA, RL, LL and C, are offered. Each cable is 1 m (3.3 ft) long. One end has a 2 mm plug for connecting to the BS4 50-9513 Isolated Preamplifier, see page I62. The other end has a 4 mm plug that connects directly to Plate and Suction Electrodes or to a Snap Fastener Connector for use with Disposable Electrodes. These Electrodes, Cables and accessories are available individually or as a complete set BS4 50-9844.

- For clinical ECG measurement for human use only
- Use only with the BS4 50-9513 Isolated Preamplifier, see page I62
- Four surface electrodes are offered:
 - A range of disposable biopotential electrodes are offered, see pages I73 to I74
 - Plate electrodes with straps
 - Suction electrode

Catalog No.	\$	Product
BS4 50-9844		5-Lead Clinical ECG Electrode Set. Complete Set of all items whose catalog numbers are shaded
BS4 50-9562		ECG Limb Cables, color-coded, 1 m (3.3 ft) leads, 2 mm plug to 4 mm plug connectors, set of 5
BS4 50-9570		4 mm Socket to Snap Fastener (press stud) Connectors, set of 5
BS4 50-5297		Plate Electrode and Strap, pkg of 1
BS4 50-9539		Plate Electrodes and Straps, set of 4
BS4 50-9604		Suction Electrode
BS4 72-7095		Disposable Foam Electrodes for Active Subjects, 38 mm (1-1/2 in), pkg of 50, (for additional Disposable Electrodes, see pages I73 to I74)



GSR Palm Input Cables and Electrodes

Both Input Cables are 1 m (3.3 ft) long and have a 3-pin DIN connector for directly connecting to the BS4 50-9505 Isolated GSR Amplifier, see page I147. One BS4 50-9521 2-Lead Input Cable with Finger Electrodes is supplied with this GSR Amplifier. These GSR Electrodes, Cables and accessories are available individually or as complete GSR Palm Disposable Electrode Set BS4 50-9810.

- For recording GSR from palm of hand and from other parts of body
- Two input cables are available:
 - 2-lead cable with attached finger electrodes
 - 2-lead cable with two 4 mm plugs to attach to BS4 50-9570 Snap Fastener Connectors for use with disposable electrodes

Catalog No.	\$	Product
BS4 50-9810		GSR Palm Disposable Electrode Set. Complete set of all items whose catalog numbers are shaded
BS4 50-9521		2-Lead Input Cable with Finger Electrodes*
BS4 50-9828		2-Lead Input Cable with 4 mm Plugs*
BS4 50-9570		4 mm Socket to Snap Fastener (Press Stud) Connectors, set of 5
BS4 72-7095		Disposable Foam Electrodes for Active Subjects, 38 mm (1-1/2 in), pkg of 50, (for additional Disposable Electrodes, see pages I73 to I74)

*Note: Pin 1 and 2 internally shorted.