



Transferpette®-12 electronic

| Capacity µl | Subdiv. µl | A* ≤ ± % µl | | CV* ≤ % µl | | Tip type | With AC adapter for | Cat. No. |
|----------------|---------------|----------------|------|---------------|------|---------------------|----------------------------------|----------|
| 0.5 - 10 | 0.01 | 1.2 | 0.12 | 0.8 | 0.08 | A, B, I***, J | Europe (continental) (230V/50Hz) | 7054 49 |
| | | | | | | | UK/Ireland (230V/50Hz) | 7054 59 |
| | | | | | | | USA/Japan (110V/50-60Hz) | 7054 69 |
| | | | | | | | Australia (240V/50Hz) | 7054 79 |
| 1 - 20 | 0.02 | 1.0 | 0.2 | 0.5 | 0.1 | A, B, J*** | Europe (continental) (230V/50Hz) | 7054 50 |
| | | | | | | | UK/Ireland (230V/50Hz) | 7054 60 |
| | | | | | | | USA/Japan (110V/50-60Hz) | 7054 70 |
| | | | | | | | Australia (240V/50Hz) | 7054 80 |
| 5 - 100 | 0.1 | 0.8 | 0.8 | 0.25 | 0.25 | C, D, K***, L, M | Europe (continental) (230V/50Hz) | 7054 53 |
| | | | | | | | UK/Ireland (230V/50Hz) | 7054 63 |
| | | | | | | | USA/Japan (110V/50-60Hz) | 7054 73 |
| | | | | | | | Australia (240V/50Hz) | 7054 83 |
| 10 - 200 | 0.2 | 0.8 | 1.6 | 0.25 | 0.5 | C, D, K***, L***, M | Europe (continental) (230V/50Hz) | 7054 54 |
| | | | | | | | UK/Ireland (230V/50Hz) | 7054 64 |
| | | | | | | | USA/Japan (110V/50-60Hz) | 7054 74 |
| | | | | | | | Australia (240V/50Hz) | 7054 84 |
| 15 - 300 | 0.5 | 0.6 | 1.8 | 0.25 | 0.75 | C***, D, L***, M*** | Europe (continental) (230V/50Hz) | 7054 56 |
| | | | | | | | UK/Ireland (230V/50Hz) | 7054 66 |
| | | | | | | | USA/Japan (110V/50-60Hz) | 7054 76 |
| | | | | | | | Australia (240V/50Hz) | 7054 86 |



* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth, steady operation. The error limits are well within the limits of DIN EN ISO 8655-2. Conformity certified to DIN 12600. A = Accuracy, CV = Coefficient of variation

** See page 56 for the definition of the tip types.

*** Tip volume less than pipettes nominal volume

Accessories

(Other accessories and spare parts can be found in the operating manual.)



Reagent reservoir

PP, high clarity.
Capacity 60 ml.
Autoclavable (121 °C).

Non-sterile, with lid.
Pack of 10.

Cat. No. 7034 59

Sterile, without lid.
Packed individually.
Pack of 100.

Cat. No. 7034 11

Sterile, without lid.
5 per bag, pack of 200.

Cat. No. 7034 09