



Liquid Handling



BRAND also offers an on-site **calibration service** (for more information, please see page 291).

**Note:**

When ordering instruments with DKD certificates, the prefix 'DKD' must be added to the order number, e.g., DKD 7047 28.

## Ordering Data

**Items supplied:**

Each Transferpette® S Digital adjustable / Fixed volume pipette is conformity certified and supplied with performance certificate and rack mount.

### Transferpette® S, Digital adjustable

Capacity µl	Description	A* ≤ ± %    µl	CV* ≤ %    µl	Subdivision µl	Tip type**	Cat. No.
0.1 - 1	D-1	2    0.02	1.2    0.012	0.001	A, I	7047 68
0.5 - 10	D-10	1    0.1	0.5    0.05	0.01	A, B, I**, J	7047 70
2 - 20	D-20****	0.8    0.16	0.4    0.08	0.02	C, D, K, L	7047 72
10 - 100	D-100	0.6    0.6	0.2    0.2	0.1	C, D, K**, L	7047 74
20 - 200	D-200****	0.6    1.2	0.2    0.4	0.2	C, D, K**, L**, M	7047 78
100 - 1000	D-1000	0.6    6	0.2    2	1	E, F, N	7047 80
500 - 5000	D-5000	0.6    30	0.2    10	5	G	7047 82
1000 - 10000	D-10000	0.6    60	0.2    20	10	H	7047 84

### Transferpette® S, Fixed volume

Capacity µl	Description	A* ≤ ± %    µl	CV* ≤ %    µl	Tip type**	Cat. No.
10	F-10	1    0.1	0.5    0.05	A, B, J	7047 08
20	F-20****	0.8    0.16	0.4    0.08	C, D, K, L	7047 16
25	F-25	0.8    0.2	0.4    0.1	C, D, L	7047 20
50	F-50	0.8    0.4	0.4    0.2	C, D, L	7047 28
100	F-100	0.6    0.6	0.2    0.2	C, D, L	7047 38
200	F-200****	0.6    1.2	0.2    0.4	C, D, M	7047 44
500	F-500	0.6    3	0.2    1	E, F, N	7047 54
1000	F-1000	0.6    6	0.2    2	E, F, N	7047 62

\* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth, steady operation. The error limits are well within the limits of DIN EN ISO 8655-2. Conformity certified to DIN 12600. A = Accuracy, CV = Coefficient of variation

\*\* See page 56 for the definition of the tip types.

\*\*\* Tip volume less than pipettes nominal volume

\*\*\*\* For use only with 2-200 µl pipette tips

## Accessories

(Other accessories for Transferpette® S pipettes can be found on page 40 and in the operating manual.)

### Starter Kit

Items supplied: 3 Transferpette® S Digital adjustable pipettes, 3 x Tip-Box (filled), 3 rack mounts.



Type	The kit includes the following Transferpette® S models	Cat. No.
MICRO	D-1, D-10, D-100	7047 90
MIDI	D-20, D-200, D-1000	7047 91
MACRO	D-1000, D-5000, D-10000	7047 92
STANDARD	D-10, D-100, D-1000	7047 93