

Ordering Data

Transferpette®, Digital adjustable

Items supplied:

Each Transferpette® Digital adjustable is conformity certified and supplied with performance certificate.

Capacity µl	A* ≤ ± % µl	CV* ≤ % µl	Subdivision µl	Cat. No.
0.1 - 1	2 0.02	1.2 0.012	0.005	7041 01
0.5 - 10	1 0.1	0.8 0.08	0.05	7041 02
2 - 20	0.8 0.16	0.4 0.08	0.1	7041 03
2 - 20	0.8 0.16	0.4 0.08	0.1	7041 04
5 - 50	0.8 0.4	0.4 0.2	0.1	7041 72
10 - 100	0.6 0.6	0.2 0.2	0.1	7041 74
20 - 200	0.6 1.2	0.2 0.4	1	7041 78
25 - 250	0.6 1.5	0.2 0.5	1	7041 76
100 - 1000	0.6 6	0.2 2	1	7041 80
500 - 5000	0.6 30	0.2 10	10	7041 82

* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth operation. A = Accuracy, CV = Coefficient of variation

Transferpette®, Fixed volume

Items supplied:

Each Transferpette® Fixed volume is conformity certified and supplied with performance certificate and calibration key.

Capacity µl	A* ≤ ± % µl	CV* ≤ % µl	Cat. No.
5	1 0.05	0.8 0.04	7041 06
10	1 0.1	0.8 0.08	7041 08
20	0.8 0.16	0.4 0.08	7041 16
25	0.8 0.2	0.4 0.1	7041 20
50	0.8 0.4	0.4 0.2	7041 28
100	0.6 0.6	0.2 0.2	7041 38
200	0.6 1.2	0.2 0.4	7041 44
200	0.6 1.2	0.2 0.4	7041 46
250	0.6 1.5	0.2 0.5	7041 48
500	0.6 3	0.2 1	7041 54
1000	0.6 6	0.2 2	7041 62
2000	0.6 12	0.2 4	7041 64

* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth operation. A = Accuracy, CV = Coefficient of variation



BRAND also offers an on-site **calibration service** (for more information, please see page 273).