


A Closer Look...

The effortless ease of operation and its unique ergonomic shape make working with the Transferpette®-8/-12 piston-operated pipette so comfortable. The use of high-quality materials makes it exceptionally light in weight, while special FKM V-rings and the stepped shape of the ejector significantly reduce the ejection force needed. Your hand remains relaxed and loose, even during a prolonged pipetting operation.



Ordering Data

Items supplied:

Each Transferpette®-8/-12 pipette is conformity certified and supplied with performance certificate, 1 Tip-Box, filled with PLASTIBRAND® pipette tips, 1 Tip-Rack for refilling, 1 stand, 1 reagent-reservoir, 1 set of sealing rings made of FKM.

Transferpette®-8

Capacity µl	A* ≤ ± %	µl	CV* ≤ %	µl	Subdivision µl	Tip type**	Cat. No.
0.5 - 10	1.6	0.16	1.0	0.1	0.05	A, B, I***, J	7036 00
2 - 20	1.0	0.2	0.6	0.12	0.1	A, B, J***	7036 02
2.5 - 25	1.0	0.25	0.6	0.15	0.1	C, D, K***, L, M	7036 04
5 - 50	0.8	0.4	0.4	0.2	0.1	C, D, K***, L, M	7036 06
10 - 100	0.8	0.8	0.3	0.3	0.1	C, D, K***, L, M	7036 08
20 - 200	0.8	1.6	0.3	0.6	1	C, D, K***, L***, M	7036 10
30 - 300	0.6	1.8	0.3	0.9	1	C***, D, L***, M***	7036 12



Transferpette®-12

Capacity µl	A* ≤ ± %	µl	CV* ≤ %	µl	Subdivision µl	Tip type**	Cat. No.
0.5 - 10	1.6	0.16	1.0	0.1	0.05	A, B, I***, J	7036 20
2 - 20	1.0	0.2	0.6	0.12	0.1	A, B, J***	7036 22
2.5 - 25	1.0	0.25	0.6	0.15	0.1	C, D, K***, L, M	7036 24
5 - 50	0.8	0.4	0.4	0.2	0.1	C, D, K***, L, M	7036 26
10 - 100	0.8	0.8	0.3	0.3	0.1	C, D, K***, L, M	7036 28
20 - 200	0.8	1.6	0.3	0.6	1	C, D, K***, L***, M	7036 30
30 - 300	0.6	1.8	0.3	0.9	1	C***, D, L***, M***	7036 32




* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth, steady operation. The error limits are within the limits of DIN EN ISO 8655-2. Conformity certified to DIN 12600. A = Accuracy, CV = Coefficient of variation

** See page 56 for the definition of the tip types.

*** Tip volume less than pipettes nominal volume

Accessories

(Other accessories and spare parts can be found in the operating manual.)



Individual stand

For 1 Transferpette®-8/-12 or 1 Transferpette® S -8/-12 pipette. Pack of 1.

Cat. No. 7034 40



Reagent reservoir,
PP, non-sterile or sterile,
please see page 55.