



9985-1005219

CPR Baby Anne

The Baby Anne manikin meets your need for a low cost, life like infant CPR trainer. It is the perfect supplement to the Resusci Baby manikin.

- Lightweight
- Affordable price for increased hands-on classroom practice
- Realistic chest compression and rise
- Inexpensive and easy-to-maintain airway system (one airway per student)
- Unique foreign body airway obstruction feature
- Head tilt/chin lift is required to open the airway

Single pack includes 1 manikin, soft pack, 6 airways, 10 foreign body practice objects, directions for use.

50x15x20 cm; 1.8 kg

E/D/F/S/J/NO/SE/I/NL

9985-1005219

CPR Baby Anne, Dark Skin

(not shown)

All other features same as 9985-1005219.

9985-1005220

CPR Baby Anne 4-Pack

All other features same as 9985-1005219. 4-Pack includes 4 manikins, soft pack, 24 airways, 40 foreign body practice objects, directions for use "CPR Baby Anne 4-Pack".

63x28x30 cm; 6.3 kg

9985-1005221

CPR Baby Anne 4-Pack, Dark Skin

4-Pack includes 4 manikins, soft pack, 24 airways, 40 foreign body practice objects, directions for use "CPR Baby Anne, 4-Pack".

63x28x30 cm; 6.3 kg

9985-1005222



9985-1005221

Replacements for 9985-1005219, 9985-1005220, 9985-1005221, 9985-1005222

Baby Anne Faces,
(pkg. 6), 25.5x21.5x12 cm;
0.79 kg

9985-1005223

Baby Anne Faces,
Dark Skin,
(pkg. 6), 25.5x21.5x12 cm;
0.79 kg

9985-1005224

Airways,
(pkg. 24), 19x16.5x15 cm;
0.35 kg

9985-1005225

Foreign Objects
8x6x2 cm; 0.01 kg

9985-1005226

**Face Connection
for Lung,**
(pkg. 50), 24x16x16 cm;
0.28 kg

9985-1005227

Infant Choking Manikin

This 9 month old baby brings to you the challenges associated with foreign body obstruction in the young paediatric age group. By placing the included obstruction in the pharyngeal cavity, you will occlude the throat and initialize the need for protocol manoeuvres to remove the object. When proper hand placement occurs with sufficient force, the object expels from the mouth. Includes two objects, clothes and carrying bag.

71x20x20 cm; 3 kg

9985-1005778



9985-1005778