



K R E A T I V E E L E K T R O P H O R E S E - T E C H N I K

PHERO-sub 1515-E



- Three tray options
- Run up to 200 samples
- Low buffer volumes
- Multichannel pipette compatible combs for speed loading

The PHERO-sub 1515-E offers a wide degree of versatility. Three tray options are available – 15 x 7 cm, 15 x 10 cm, and 15 x 15 cm – allowing the choice of one, two or all

three gel length options at the time of purchase. Further purchases of additional accessories are no longer required. Maximising comb and tray options allow up to 210 samples to be resolved per gel. The 15 cm total run length allows restriction fragment or other close MW sample bands to be easily separated and identified. Speed loading is accomplished using 10, 14, 18, 28 and 30 sample multichannel pipette compatible combs.

Technical specifications

<i>Gel dimensions (W x L)</i>	15 x 7 cm 15 x 10 cm 15 x 15 cm
<i>Unit dimensions</i>	26.5 x 17.5 x 9 cm
<i>Max. Sample capacity</i>	15 x 7 cm tray – 70 samples 15 x 10 cm tray – 140 samples 15 x 15 cm tray 210 samples
<i>Buffer volume</i>	500 ml
<i>Combs</i>	1, 2, 4, 10, 10(MC), 12, 14(MC), 16, 18(MC), 20, 28(MC), 30(MC), 35
<i>Thicknesses</i>	0.75, 1, 1.5, 2 mm



Loading guides