


PEAK EXPIRATORY FLOW RATE - NORMAL VALUES

For use with EU/EN13826 scale PEF meters only



In men, readings up to 100 L/min lower than predicted are within normal limits. For women, the equivalent figure is 85 L/min. Values are derived from Caucasian populations.

Mini-Wright
(Standard Range)
EU scale
(EN 13826)

Blue text
on a yellow
background



Single Patient Use
Part Ref: 3103388

Multiple Patient Use
Part Ref: 3103387

NHS Logistics
Code : FDD 609

Mini-Wright
(Low Range)
EU scale

Blue text
on a yellow
background



Single Patient Use
Part Ref: 3104708

Multiple Patient Use
Part Ref: 3104710

Clement Clarke has developed mathematical equations that will allow conversion of P.E.F. readings from Wright-McKerrow scale to EN 13826 scale, and vice-versa. Contact us directly, or visit the website



Connect With Us

