



LIFE/FORM® TRACHEOSTOMY CARE SIMULATOR

- Suctioning
- Cuff Inflation
- Dressing Changes
- Visual Landmarks
- Visual Verification of Bronchial Location

The Life/form® Tracheostomy Care Simulator is the most realistic approximation of a patient available today. Careful attention to detail is incorporated into this pseudo-patient. Suctioning techniques, proper cuff inflation, dressing changes and other techniques can be practiced and critiqued. The simulator replicates an adult male containing oral and nasal passages and all appropriate anatomy - pharynx, epiglottis, trachea, esophagus, stoma, cricoid cartilage and representative cervical vertebrae. The chest area contains left and right bronchi and the bronchial tree. The esophagus dead ends 2" below the trachea opening. The side of the neck and chest area over the bronchial tree allow viewing of the suction and trachea tube. The viewing window at the bottom of the simulator shows the location of a suction catheter in the right bronchus (seen) or left bronchus (unseen). Oral, nasopharyngeal, nasotracheal and tracheal suctioning can be practiced. Cleansing of the stoma area can be performed just as on an actual patient. Changing dressings and ties is also possible and the experience is very close to an actual one. Confidence and practice are at your fingertips with this Life/form® Tracheostomy Care Simulator. Comes complete with instruction guide, hard carrying case and REN cleaner. Tracheostomy tube is not included.

CODE	SIZE	SHIPPING WEIGHT
LIFE/FORM® TRACHEOSTOMY CARE SIMULATOR		
LF01083	53 x 53 x 28 cm	11.3kg

LIFE/FORM® PATIENT EDUCATION TRACHEOSTOMY CARE SET

- Adult and Infant Manikins
- Suctioning
- Dressing Changes
- Cuff Inflation

The Life/form® Patient Education Tracheostomy Care Set is ideal for teaching patients and caregivers the skills they will need to perform at home. The Adult Manikin includes oral and nasal passages, the pharynx, epiglottis, trachea, esophagus, tracheostomy site and cricoid cartilage. The mouth and jaw are flexible to allow oral suctioning. The trachea branches into the left and right bronchial trees. Students can practice suctioning techniques, proper cuff inflation, dressing changes and other techniques. Cleansing and maintenance of the external tracheal site can be performed just as on an actual patient. The Infant Manikin can be used to teach parents and caregivers basic tracheostomy care skills and is used extensively with children to comfort them and improve their self-esteem. The Life/form® Tracheostomy Care simulator set comes complete with one adult and one infant manikin, one adult and one infant tracheostomy tube, instruction guide and hard carrying case.

CODE	SIZE	SHIPPING WEIGHT
LIFE/FORM® PATIENT EDUCATION TRACHEOSTOMY CARE SET		
LF01159	-	12.2kg
ADULT TRACHEOSTOMY CARE MANIKIN		
Adult manikin comes with tracheostomy tube, instruction guide and hard carrying case. One-year warranty		
LF01168	-	6.8kg
INFANT TRACHEOSTOMY CARE MANIKIN		
Infant manikin comes with tracheostomy tube, instruction guide and hard carrying case. One-year warranty		
LF01167	-	5.9kg



Connect With Us

