



Rapid Test Formats

Multi-Drug Oral Fluid Test Cassette

Features:

A rapid one step test for the qualitative detection of the following drugs of abuse and their principal metabolites in human oral fluid at specified cut off level:

- Amphetamine (AMP),
- Benzodiazepines (BZO),
- Cocaine (COC),
- Marijuana (THC),
- Methamphetamine (MET),
- Morphine (MOP).



Handheld DOA Reader

Quantitative result within several minutes
Compatible with Wondfo DOA products
Portable device in one hand



Cup Tests

DOA Urine Cup Test is a popular format for maximum 18 drug tests one time. Customers can combine any drug test items to this cup and get the result from a single urine sample. This formats contain the incorporated temperature strip which can protect the prompt urine collection, check urine temperature. No need to handle any urine during the test. Tester can get the result in 3~5 minutes from the showing window. The test result is over 99% accuracy. This format is very easy to use and no special training is needed. The instruction insert is applied to each test.

