

## » 净化工作台

### PURIFICATION TABLE



S·SW-CJ-1F单人双面净化工作台

#### 垂直层流形标准型净化工作台

参 数	型 号	S.SW-CJ-1F S.SW-CJ-1FD	S.SW-CJ-2F S.SW-CJ-2FD	S.VS-1300 S.VS-1300U
洁净等级	100级@≥0.5um(联邦209E)			
菌 落 数	≤0.5个/皿 · 时(Φ90mm培养平皿)			
平均风速	0.3 ~ 0.6m/s(推荐使用风速0.5m/s)			
噪 音	≤62dB (A)			
振动半峰值	≤5um	≤5um	≤4um	
照 度	≥300Lx			
电 源	AC 单相220V/50Hz			
最大功耗	0.4Kw	0.8Kw	0.8Kw	
净重/毛重	80/130Kg	130/210Kg	130/210Kg	
规 格 尺 寸 3	W1:工作区宽	87	136	112
	D1:工作区深	70	70	65
	H1:工作区高	52	46	72
	W:装置外形宽	99	148	136
	D:装置外形深	72	72	80
	H:装置外形高	166	166	180
包装箱尺寸cm	114 × 600 × 186	164 × 86 × 168	144 × 94 × 191	
高效过滤器规格、数量	820 × 600 × 60 × ①	610 × 610 × 69 × ②	610 × 610 × 69 × ②	
萤光灯/紫外灯规格、数量	20W × ①/20W × ①	20W × ①/20W × ②	20W × ①/20W × ②	
适用人数	1F:单人双面 1FD:单人单面	2F: 双人双面 2FD:双人单面	双人单面	
注	U为带紫外灯			

#### 适用范围:

广泛适用于医疗卫生、制药、化学实验、电子、国防、精密仪器、仪表等行业。

#### 结构特点:

- 垂直流形,准闭合式台面,可有效防止外部气流诱入及操作异味对人的刺激。
- S.VS-1300配以弹簧式上下移门,定位灵活方便。
- 采用进口不锈钢台面和后衬板,配以玻璃侧档板,使整个工作区显得宽敞明高。
- 铝合金制成的风门可任意调节预过滤进风风量。
- 采用可调风量送系统,轻触型开关调节风速,使工作区风速处于理想状态。

#### USAGE:

The product is widely used for medical hygiene pharmacy, chemical test, electricity, national defence precision instruments and apparatus.

#### STRUCTURAL FEATURES:

- A. Horizontal single-direction flow
  - It's horizontal flow and open-top, to operate conveniently.
  - The towedeh switching may regnlate the power. The wind machinery system regulated volume guarantees the working room in ideal state.
- B. Vertical single-direction flow
  - It's vertical flow and para-closed top, to prevent the body from stimulation of the strang operation smell, and from the ontside gas.



S·VS-1300型双人单面垂直净化工作台