

## **Abdominal Open and Closure Trainer**

The Advance Abdominal Open and Closure Skin Pad has been developed to simulate real skin tissue and for ease of use by course organisers doing basic surgical courses. It is supplied with a pink skin surface with subdermis for intra-dermal suturing. A white Linear Alba makes up a further layer.

- Easy maintenance, requiring only a set of balloons and set of skins
- Up to 4 good incisions can be made from each set of skins
- · Skins cut and suture well

W19356

W19359

W19360

Carrying Bag, 17x36x30 cm; 4.8 kg

- Included with the kit are four balloons and a set of skins (1 red, 1 white)
- · Internal planes within the unit prevent the balloon from moving around
- Rubber suction feet secure the base providing easy working conditions

W19330

**Replacements for W19330** 

4 Balloons and 2 Skins

W19331





## W19331



FOR ADVANCED TRAUMA LIFESUPPORT COURSES

This innovative design not only simulates realistic body tissue but provides an alternative to the use of animals. The design is unique in having a replaceable chest wall with ribs that are encased into the realistic simulated flesh that forms the torso. This fleshy part is replaceable and can used for up to 25 incisions. It is housed in a white plastic torso-shaped box that is firmly secured and has strong sucker-feet to ensure non slip use on smooth table-top surfaces. The torso presents as a patient lying on his back.

- · Surgical access is lateral to the pectoralis major
- The arm has been removed for easier access and extended for anatomical
- The ribs can be clearly felt below the surface of the skin
- Trainees can make an incision through the fleshy part, surgically dividing the tissue with blunt forceps until the plural cavity is clearly felt by finger
- Complete finger rotation is possible, allowing the trainee to ensure that there are no obstructions before a drainage system can be introduced
- The popping effect of passing through the pleura is realistic and a unique feature
- A drain can be securely sutured into position onto this simulated flesh
- The trainer needs no preparation, no refrigeration, no messy disposal after use and no unpleasant odours

Trainer includes a replaceable:

- · Chest wall
- Pericardiocentesis where fluid can be drawn into a syringe
- A realistic Tension Pneumothorax where air can clearly be heard and felt exhaling between the 2nd and 3rd ribs

57x37x43 cm; 4.8 kg

□ E



Connect With Us



















W19359

