

- ✓ XC-406-1 / 406A-1 Whole Body Basic CPR Manikin Style 100 (Male / Female)
- ✓ XC-406-2 / 406A-2 Whole Body Basic CPR Manikin Style 200 (Male / Female)
- ✓ XC-406-5 / 406A-5 Whole Body Basic CPR Manikin Style 500 (Male / Female)



XC-406-1

XC-406-2

XC-406-5

Implement standard: AHA (American Heart Association) guideline for CPR

Feature:

CPR manikin feature:

- 1) Clearly anatomic characteristic, realistic touch feeling and lifelike skin color, vivid appearance.
- 2) Simulate vital signs.
 - (a) Pupil state: contrast observation of pupils, one is dilated, the other one is contracted.
 - (b) Carotid artery response: simulate carotid artery pulse.
- 3) Open airway.
- 4) Can perform artificial breathing and external cardiac compression.

CPR instrument feature:

- 1) Available CPR training.
- 2) Electronic monitoring of compression site.
- 3) Indicators show inflation volume, correct inflation volume: 800-1200ml.
 - (a) Insufficient inflation volume, yellow indicator.
 - (b) Proper inflation volume, green indicator.
 - (c) Excessive inflation volume, red indicator.
- 4) Indicators show compression depth: correct compression depth: 4-5cm.
 - (a) Insufficient compression depth, yellow indicator.
 - (b) Proper compression depth, green indicator.
- 5) Voice Prompt.
- 6) Operation rate: ≥ 100 times/min.
- 7) Operation cycle: 2 valid inflation after 30 valid compression, 5 cycles.
- 8) Set the operating time, values in seconds.
- 9) Result Printing: XC-406-2/XC-406-5.

Packing: 1pcs/carton, 95x34x58cm, 21kgs

