

Anesthesia Masks and Inhalation Narcosis Chamber



Previous



Next

Table of Contents

Product Index

Search

Home

Contact Us



Previous



Next



Cat Ventilation/Anesthesia Masks

These soft pliable vinyl plastic masks conform to the facial contours of a cat. They are fitted with a standard 15 mm (3/5 in) taper for connection to anesthesia machines, see pages F16 to F19. They are available in three sizes and are supplied individually or as a convenient set containing one mask of each size.

The mask diameters are 25 mm (1 in), 38 mm (1-1/2 in), and 50 mm (2 in) for the small, medium and large size masks respectively.

Catalog No.	\$	Product
BS4 59-8201		Cat Ventilation/Anesthesia Mask, Small
BS4 59-8202		Cat Ventilation/Anesthesia Mask, Medium
BS4 59-8203		Cat Ventilation/Anesthesia Mask, Large
BS4 59-8204		Set of 3 Masks, one of each size



Inhalation Narcosis Chamber

A clear acrylic hinged cover permits observation of animals. All other components are made of polypropylene and are not affected by anesthetic compounds.

A folding perforated platform is removed to fill the bottom with cotton and replaced to support the animal. Anesthetic is poured into the small funnel at the side and runs into the bottom of the chamber where it is absorbed by cotton. A second stopcock on the other side of the chamber permits the entrance of fresh air during anesthesia to prevent suffocation.

Maximum ID is 197 mm (7-3/4 in); top ID is 178 mm (7 in); height from platform to top is 127 mm (5 in); height under platform is 57 mm (2-1/4 in); overall height is 197 mm (7-3/4 in).

Catalog No.	\$	Product
BS4 59-6717		Inhalation Narcosis Chamber

Connect With Us



Canine Ventilation/Anesthesia Masks

These masks provide a convenient way of providing artificial ventilation or administration of anesthetic gases to dogs. They are constructed of rugged transparent polycarbonate. One end has a standard 15 mm (3/5 in) taper for connection to anesthesia tubing, while the other end has a soft flexible rubber diaphragm that molds itself to the contours of a dog's snout.

Two sizes are available. The small mask has an unstretched diaphragm opening of 2.5 cm (1 in), while the large mask has an unstretched opening of 4.4 cm (1-3/4 in).

Catalog No.	\$	Product
BS4 59-7976		Canine Ventilation/Anesthesia Mask, Small
BS4 59-7977		Canine Ventilation/Anesthesia Mask, Large
BS4 59-7978		Replacement Diaphragm for Small Mask
BS4 59-7979		Replacement Diaphragm for Large Mask



- For rats, guinea pigs, birds, and mice
- Simple, yet effective design

Small Animal Anesthesia Masks

These Anesthesia Masks were designed with smaller animals in mind. They have a simple cylindrical shape with a rubber diaphragm. A small slit is cut in the diaphragm to accommodate the animal's snout or beak. The opposite end tapers down to a standard press fitting port measuring 15.3 mm OD x 11.9 mm ID. These Masks can be used with bain circuits, non-rebreathing systems, or 'Y' pieces. Includes 1 diaphragm.

Catalog No.	\$	Product
BS4 59-8254		Mouse Anesthesia Mask
BS4 59-8255		Rat, Guinea Pig and Bird Anesthesia Mask
BS4 72-3923		Replacement Diaphragms, pkg. of 10