

**Condom Training Model, White Skin Tone**

This model of an erect penis with testicles is used to quickly learn how to use a condom safely. The anatomical structures and its firmness are absolutely realistic, so that your students can train putting on and removing a condom in a realistic way. Supplied with 12 dry training condoms and a carrying bag. 7.5x7.5x19.5 cm; 0.35 kg

E/D/S/F/P/I/J www.3b.com

L42

Condom Training Model, Coloured Skin Tone

Like L42.

L43

**Condom Training Model**

Demonstrate the proper use of condoms by using this realistic model. Consists of an erect penis, 12 condoms, syringe and artificial semen (UV-fluorescent fluid) to simulate ejaculation. Mounted on a stand with suction cups and supplied with carrying bag. 35.5x15x16.5 cm; 2.3 kg

E

W43001

**Condom Training Models**

These low priced Styrofoam models are excellent for use in larger group lessons and can be reused several times. The models can be fixed to the desktop with adhesive tape, so that both hands are free for rolling the condom into position. Supplied without condoms. 14.5 cm

W19101

**Training Model for a Female Condom**

This model shows the labia and vagina up to the cervix in a simplified representation for didactic reasons, and is used for demonstrating and learning the insertion of a female condom. Supplied with one Femidom condom. 12 cm; 0.15 kg

L41/2

Dark skin

L41/1

**AIDS Virus**

This model of the HI-Virus, enlarged millions of times, shows the outer lipid membrane with protein structures, and the internal nucleus which contains the viral hereditary matter (RNA). The nucleus is removable and condoms can be put inside to provide a message regarding measures to take in protecting against HIV. Supplied without condoms. Mounted on base. 18x13x13 cm; 0.7 kg

L40

Family Planning Educator

Desktop simulator for training and demonstration of: introduction and removal of a diaphragm, an IUD or sponge contraception devices; normal and abnormal uterine positions; bi-manual examination technique. Supplied with: one anteverted uterus with clear upper half to illustrate correct position of IUD; one uterus to illustrate normal anteversion and retroversion; cervix with patent os attaches to uterus suspended within pelvic cavity; soft plastic stomach cover; carrying bag.

E

W45009

**Options and Replacements for W45009 and W43001****Artificial Semen (UV-fluorescent fluid)**

W43002

Options and Replacements for L42, L43 and W19101**12 Dry Condoms**

W43003

Options and Replacements for W43001**12 Dry Condoms**

W43003

I.U.D Trainer

Hand held trainer which is a suitable aid for understanding correct positioning of I.U.D. (Intrauterine Device) in the uterus. Made of durable plastic, the trainer features a transparent cover which allows easy visualization of insertion and placement of I.U.D. (I.U.D. not included). 6x40x45 cm

E

W44615



W45154

Female Condom Model

This model represents an anteverted uterus in a simulated pelvic cavity with soft vulva and vagina and a soft plastic stomach cover. A perfect tool to demonstrate the use of a female condom, contraceptive sponge and cervical cap. Supplied with carrying bag. 17.8x24x14 cm

E

W45154

I.U.D Trainer

This anatomically accurate model represents a section of the uterus, ovaries and fimbriae. The uterus is covered by a clear plastic window to allow easy sight of insertion and placement of I.U.D. (I.U.D. not included).

E

W45152



W45152