



Hitachi 902

Reliability, quality and convenience make the Roche/Hitachi 902 the clinical chemistry workhorse for the small workload laboratory.

Compact, sample selective clinical chemistry analyser with a capacity up to 36 tests on board (including 3 ISE) and a throughput up to 200 tests per hour respectively up to 300 tests with optional ISE

A compact clinical chemistry analyzer for laboratories with small workloads.

The Roche/Hitachi 902 has an on-board capacity of 26 tests (including 3 ISEs – Na⁺, K⁺, Cl⁻), throughput of up to 200 tests per hour and up to 300 tests with optional ISE unit.

For analysis of serum, plasma, urine, cerebrospinal fluid (CSF) and hemolysate.

Roche/Hitachi 912

Analyzer for laboratories with small to medium workloads.

The Roche/Hitachi 912 is a flexible system for both classical and special clinical chemistries, homogeneous immunoassays including STAT analyses. It has an on-board capacity of 35 tests (including 3 ISEs – Na⁺, K⁺, Cl⁻) and throughput of up to 360 tests per hour (up to 720 tests with optional ISE unit).

For analysis of serum, plasma, urine, cerebrospinal fluid (CSF) and hemolysate.

Customer benefits

- **Broad menu for test consolidation**
The Roche/Hitachi test menu includes assays for routine clinical chemistry, protein profiling, therapeutic drug monitoring and drugs of abuse screening (Roche/Hitachi 912 only). With up to 26 (Roche/Hitachi 902) or 35 (Roche/Hitachi 912) on-board assays, the systems can consolidate routine testing or special testing in a variety of laboratory environments.
- **Proven reliability and quality**
Roche/Hitachi systems have set the standard in the diagnostic industry. Proven and dependable Roche/Hitachi systems ensure high-quality results 24 hours a day in laboratories worldwide.

- **Efficient operation and workflow**
Simplified operation combined with maximum flexibility and security is achieved by such features as user-friendly touch-screen and bar-coded sample identification.

Additional feature

- **Improved lab efficiency**
With auto-calibration, automatic pre-dilution and rerun, with the Roche/Hitachi 912 system the laboratory workload can be processed with minimal user intervention.