Induction Chamber, Anesthetizing Boxes and **Backdraft Anesthesia Table**



- Clear acrylic for easy animal observation
- Two sizes
- · Ideal for small animals
- · Stainless steel false bottom for easy cleaning

Induction Chamber

This cylindrical acrylate chamber is ideal for anesthetizing small animals such as rodents, guinea

pigs, and birds. It has two, 1/4 inch barbed hose fittings. One for gas delivery and one for gas scavenging. The gas delivery barb directs the flow of anesthetic gas towards the bottom of the chamber. The chamber has a routed top that ensures a perfect fit. There are no corners so cleaning is easy. Includes 10 ft. of 0.25" ID tubing.

_									
D	ı	m	e	n	S	ı	O	n	S

Inner, H x Dia	20.3 x 20.3 cm (8 x 8 in)	
Small	8.89 x 13.3 cm (3.5 x 5.25 in)	
Large	20.3 x 20.3 cm (8 x 8 in)	
Overall, H x Dia	30.5 x 27.9 cm (12 x 11 in)	
Small	17.8 x 17.8 cm (7 x 7 in)	
Large	30.5 x 27.9 cm (12 x 11 in)	
Barbs	Accepts, ID x OD, 6 x 8 mm (1/8 x 1/4 in) tubing	

Order# **Product**

DC1 /2-/110	induction Chamber, Small
DC1 60-5246	Induction Chamber, Large



Anesthetizing Boxes

These anesthetizing boxes are constructed of 6 mm (1/4 in) thick, transparent plastic, with all joints securely cemented. The transparent plastic allows the animal to be under constant observation. Two 4.8 mm (3/16 in) nylon tubing connectors are

cemented on opposite ends of the box, one at the top and the other at the bottom, for introducing the gas.

These boxes are available in two sizes: small and large. The small box is for rodents. The top slides open lengthwise on machined grooves. The large box is for larger animals such as rabbits and cats. Its top has three hinges for easy opening. Two fasteners squeeze the top down against a rubber gasket to seal the joints.

The small box measures, H x W x D, 10.7 x 25.7 x 11 cm (4 x 10 x 4.5 in). The large box measures, H x W x D, 29 x 44 x 29 cm (11.5 x 19.5 x 11.5 in).

*Note: Boxes cannot be cleaned with alcohol or they will crack.

Order# **Product**

DC1 50-0108	Anesthetizing Box, Small
DC1 50-0116	Anesthetizing Box, Large

Waste gases are carried away from the table through a 3 inch flexible

tube which must be vented to the outside. Fifteen feet of the tubing is

provided with the table. The table rests on 4 large swivel locking casters. It features a retractable power cord and size 'E' anesthetic gas tank

holders. The on/off switch and duplex AC power outlet are located on the



- · Pulls anesthetic waste gases down and away from surgeon
- · Better access to animal than conventional fume hoods
- One-piece table top construction, prevents debris accumulation

Specifications

backsplash of the table top.

- Promoduono			
Dimensions, H x W x D:			
Overall	106.7 x 78.7 x 127 cm (42 x 31 x 50 in)		
Work Surface	91.4 x 55.9 x 111.7 cm		
Dimensions	(36 x 22 x 44 in)		
Power Supply	110 V		

Backdraft Anesthesia Table

The Backdraft Anesthesia Table is a mobile surgical station used to pull anesthetic gases away from the surgeon. The open top design and table top work height mean that the typical workspace restrictions of fume hoods are avoided. The high air flow (155 cubic feet/ minute) provides safety for the researcher. The sturdy one-piece construction of the 16 g type 304 stainless steel work top is sanitary and easy to clean.

Connect With Us















Order# **Product**

DC1 72-2625 Back Draft Anesthesia Table

