



TANDA IID

Model	Voltage/Power	ChamberSize	OutsideSize	PackingSize	NetWeight
TANDA IID18	220V/1500W	245X320mm	441*586*395mm	660*525*495mm	47.4kg
TANDA IID23	220V/1700W	245X450mm	441*696*395mm	770*530*495mm	53.4kg

- The new design model with opening water tank on the top. Strictly meet the standard of EN13060, which makes wash and cleaning much more convenient.
- It adopts European B standard with 3-times pre-vacuum preceding vacuum drying.
- Vacuum measurement reaches -0.8bar. It sterilizes and sanitizes miscellaneous packed or unpacked medical instruments. Whether solid or A-type hollow, whether with multi-apertures or inset pipes. It efficiently tidy air from every hollow area of any hollow instrument, such as hand-piece, making sure that water steam reaches every corner.
- It is installed with BOWIC&DICK which measures the penetration of water steam.
- All-round computerized fuzzy control. Digital display modularized panel for easy operation.
- Mini printer can be attached to record the process of sterilization.