

Tenbroeck Tissue Grinders, Wheaton

Supplier: WHEATON



For homogenizing delicate tissues such as intestines, liver, and heart. Designed for hand operation, grinders can also be slow-speed motor driven by inserting a rubber stopper (size 00) and metal rod into the pestle. Hollow handle permits packing with ice. Grinders also feature an expanded reservoir and pouring lip.

Interchangeable pestles and tubes are made of Wheaton 33® low extractable borosilicate glass. Clearance between pestle and grinding chamber is 0.09±0.16mm (0.0035±0.0065"). Glass meets ASTM Type 1 Class A and USP Type 1 standards.

Click to enlarge



[TISSUE GRINDER TENBRO 15ML CS2](#)

VWR Catalog Number	Unit
62400-530	Case of 2
14231-358	Case of 2
62400-493	Case of 2
62400-551	Case of 2
62400-518	Case of 2

	Capacity	Chamber Dimensions	Overall Length	Supplier No.
	15 mL	22O.D. x 94L mm	250 mm	357426
	1 mL	11O.D. x 48L mm	140 mm	357421
	2 mL	11O.D. x 50L mm	140 mm	357422
	40 mL	32O.D. x 140L mm	320 mm	357428
	7 mL	16O.D. x 82L mm	190 mm	357424



Connect With Us

