

# OPTIMA MACINTOSH

سوماتكو  
SOMATCO  
www.somatco.com



Optima™  
Timesco England




THE MACINTOSH BLADE, DEVELOPED BY R. MACINTOSH IN 1941 IS THE MOST POPULAR BLADE IN USE FOR INTUBATION OF ADULT PATIENTS UNDER NORMAL CIRCUMSTANCES.

THE CURVATURE OF THE BLADE ALLOWS THE TIP TO FALL NATURALLY INTO POSITION IN THE VALLECULA OF THE PATIENT AND THE WIDE FLANGE ASSISTS IN HOLDING THE TONGUE SAFELY ASIDE DURING INTUBATION.




2960.150.25  
LARGE ADULT No.4  
155x18MM




2960.150.20  
ADULT No.3  
130x18MM



2960.150.15  
CHILD No.2  
100x17MM




2960.150.10  
INFANT No.1  
92x14MM



2960.150.05  
NEONATE No.0  
80x14MM





**O**ptima - one brand, one standard, one choice: In our opinion providing the optimum in performance, ease of maintenance and blade variety of any laryngoscope system currently available.

Optima is our flagship brand offering a premium quality fibre optic laryngoscope system at a competitive price. It's major benefits include:

*Fibre optic illumination* - providing the user with a bright, white light (courtesy of a 2.8v halogen bulb) which is directed exactly to where illumination is needed by the fibre optic core.

*ISO fitting* - All Optima, ViewMax and Eclipse blades are made to the ISO 7376/3 standard fitting for fibre optic laryngoscopes ensuring that all blades and handles in the system are completely interchangeable.

*Modular Design* - Optima is a modular system, breaking down into just seven major components, any of which can be easily replaced in case of malfunction.

*Replaceable Fibre Optics* - The ability to quickly remove the fibre optic cores for cleaning, maintenance or replacement enables the system to be simple and cost effective to maintain.

*Versatility* - The Optima system features a range of over a dozen different blade patterns, covering specialised areas such as difficult intubation, paediatrics and cleft palate work. Disposable blades are available to compliment the range in cases where cross infection may be of particularly high risk. Also, the handles are available with the option of battery power or long lasting rechargeable cells.



OPTIMA BLADES AND HANDLES ARE COLOUR CODED WITH GREEN BLOCKS AND BANDS TO ENABLE THEM TO BE EASILY IDENTIFIED AS COMPATIBLE.



THE OPTIMA LEFT HANDED MACINTOSH BLADE HAS THE FLANGE ON THE OPPOSITE SIDE FROM THE REGULAR MACINTOSH. AS WELL AS FOR THE COMFORT OF THE LEFT HANDED USER, THIS DESIGN IS ALSO USEFUL FOR PATIENTS WITH ABNORMALITIES OF THE RIGHT SIDE OF THE FACE OR OROPHARYNX, INTUBATING IN THE RIGHT LATERAL POSITION OR POSITIONING A TUBE DIRECTLY ON THE LEFT SIDE OF THE MOUTH.



2960.170.20  
LEFT HANDED  
ADULT No.3  
130x18MM

(2960.170.25 LEFT HAND ADULT No.4  
155x18MM ALSO AVAILABLE)

Connect With Us

