Restraint Devices







Table of

Contents

Product

Index

Search

Home



- Perfect for injections or blood sampling
- Entire tail is completely exposed for injections or sampling

These popular Broome Restrainers provide fast and easy animal insertion. The cylindrical design of these Restrainers solves the problem of the rodents turning once they are introduced into the Restrainer.

To place the animal into the Restrainer, loosen white Delrin nose piece and grasp the animal by the tail. Slide the animal into the Restrainer by pulling its tail along the open groove until the animal is completely enclosed in the cylinder. Then slide the nose piece into the groove on the opposite side of the Restrainer and place it up to the nose of the animal. Acrylic. Not autoclavable.

| Broome Rodent Restrainers | | | | | | | |
|---------------------------|----|------------------------|--------------------------|---------------------|--|--|--|
| Catalog No. | \$ | Animal Weight Range | Restraining Barrel ID | Length Min x Max | | | |
| BS4 52-0460 | | 15 to 30 g | 25.4 mm | 41 x 83 mm | | | |
| BS4 52-0470 | | 30 to 70 g | 31.75 mm | 50 x 100 mm | | | |
| BS4 52-0478 | | 70 to 125 g | 38 mm | 73 x 124 mm | | | |
| BS4 52-0486 | | 125 to 250 g | 50.8 mm | 70 x 178 mm | | | |
| BS4 52-0494 | | 250 to 500 g | 63.5 mm | 58 x 195 mm | | | |

Universal Rodent Restrainers







Next

- These Universal Rodent Restrainers have a unique raised inner floor which provides a better fit for smaller rodents and which may be removed for larger rodents. The inner floor is attached to the bottom by two nylon screws. Knurled, nickel-plated brass screws
- Restrains mice down to 10 grams!
- Removable raised inner floor
- Adjustable headgate and removable tailgate
- Headgate fits into any of four positions for various sizes of rodents
- Tailgate allows maximum closeness of occluding cuff to body when testing indirect blood pressure

serve as base fasteners and legs. These Universal Rodent Restrainers are constructed of 6.4 mm (1/4 in) thick, rigid, clear acrylic plastic. Openings in the headgate reduce restrictions on the rodent's breathing. Not autoclavable.

| Universal Rodent Restrainers | | | | | | | |
|------------------------------|----|----------------------------------|-------------------------------------|-------------------------------------|--|--|--|
| Catalog No. | \$ | Animal Weight Range | Restraining Area Minimum | ı, H x L Maximum | | | |
| BS4 56-4499 | | 10 to 40 g (0.4 to 1.4 oz) | 28.6 x 51 mm (1-1/8 x 2 in) | 41.3 x 102 mm (1-5/8 x 4 in) | | | |
| BS4 56-4500 | | 150 to 300 g (5.3 to 10.6 oz) | 41.3 x 124 mm (1-5/8 x 4-7/8 in) | 54 x 171 mm (2-1/8 x 6-3/4 in) | | | |
| BS4 56-4501 | | 250 to 600 g (8.8 to 21.2 oz) | 47.6 x 152 mm (1-7/8 x 6 in) | 60.3 x 216 mm (2-3/8 x 8-1/2 in) | | | |

💷 Rodent Restraint Bags



- Variety of sizes and shapes
- Multiple colors
- Safe, humane and disposable
- Includes twist ties to secure bags
- Multi-use

These Rodent Restraint Bags constructed of Easy-Flex-PE^m provide effortless wrapping and added animal comfort. Pre-cut sizing eliminates the need for cutting down larger bags for smaller animals, and the optional Body-Contour^m Cone shape provides easier animal capture and restraint and increased animal comfort. Optional Black and Rose color bags reduce natural light to provide an environment that aids in reducing animal agitation and stress. Restraint Bags can be used in a variety of procedures including tail tattooing and injections, IP injections, oral dosing and decapitation.

| Rodent Restraint Bags | | | | | | | | | | |
|---|-------------------|------------------|----------------|-------------|--|--|--|--|--|--|
| Body-Contour [™] Shape, pkg. of 40 | | | | | | | | | | |
| Clear Cat. No. | Black Cat. No. | Rose Cat. No. | Animal Size | Bag Size | | | | | | |
| BS4 72-6400 \$ | BS4 72-6401 | BS4 72-6402 | Small | 120 g | | | | | | |
| BS4 72-6403 \$ | BS4 72-6404 | BS4 72-6405 | Medium | 225 g | | | | | | |
| BS4 72-6406 \$ | BS4 72-6407 | BS4 72-6408 | Large | 450 g | | | | | | |
| Traditional Triangle Shape, pkg. of 50 | | | | | | | | | | |
| BS4 72-6409 \$ | - | _ | X-small | < 20 g | | | | | | |
| BS4 72-6410 \$ | BS4 72-6411 | BS4 72-6412 | Small | 120 g | | | | | | |
| BS4 72-6413 \$ | BS4 72-6414 | BS4 72-6415 | Medium | 225 g | | | | | | |
| BS4 72-6416 \$ | BS4 72-6417 | BS4 72-6418 | Large | 450 g | | | | | | |



Connect With Us

