

Erlenmeyerkolben · Erlenmeyer flasks

Erlenmeyerkolben, Borosilikatglas 3.3, mit Bördelrand, graduert

Erlenmeyer flasks, borosilicate glass 3.3, with beaded rim, graduated

Cat. No.	Inhalt Capacity	Ø Boden x Ø Hals Ø bottom x Ø neck	Höhe Height	unit
41 102 02	25 ml	42 x 22 mm	75 mm	10
41 102 03	50 ml	51 x 22 mm	90 mm	10
41 102 04	100 ml	64 x 22 mm	105 mm	10
41 102 06	200 ml	79 x 34 mm	135 mm	10
41 102 07	250 ml	85 x 34 mm	145 mm	10
41 102 08	300 ml	87 x 34 mm	160 mm	10
41 102 10	500 ml	105 x 34 mm	180 mm	10
41 102 13	1000 ml	131 x 42 mm	220 mm	10
41 102 14	2000 ml	166 x 50 mm	280 mm	10
41 102 15	3000 ml	187 x 50 mm	310 mm	4
41 102 17	5000 ml	220 x 50 mm	365 mm	1

Enghals
Narrow neck
DIN ISO 1773

41 103 02	25 ml	44 x 32 mm	70 mm	10
41 103 03	50 ml	51 x 34 mm	85 mm	10
41 103 04	100 ml	64 x 34 mm	110 mm	10
41 103 06	200 ml	79 x 50 mm	131 mm	10
41 103 07	250 ml	85 x 50 mm	140 mm	10
41 103 08	300 ml	87 x 50 mm	156 mm	10
41 103 10	500 ml	105 x 50 mm	175 mm	10
41 103 13	1000 ml	131 x 50 mm	220 mm	10
41 103 14	2000 ml	165 x 70 mm	280 mm	10

Weithals
Wide neck
DIN 12 385

Erlenmeyerkolben, Borosilikatglas 3.3, mit Normschliff, graduert

Erlenmeyer flasks, borosilicate glass 3.3, with standard ground joint, graduated

Cat. No.	Inhalt Capacity	Normschliff Standard ground joint	unit
43 103 02	25 ml	NS 14/23	10
43 104 03	50 ml	NS 19/26	10
43 105 04	100 ml	NS 24/29	10
43 106 04	100 ml	NS 29/32	10
43 105 07	250 ml	NS 24/29	10
43 106 07	250 ml	NS 29/32	10
43 106 10	500 ml	NS 29/32	10
43 106 13	1000 ml	NS 29/32	10

DIN 12 387

Normschliff-Stopfen, Borosilikatglas Duran, hohl, mit Sechskantgriff und geschlossener Spitze

Standard ground-joint stoppers, borosilicate glass Duran, hollow, with hexagonal grip and closed tip

Cat. No.	unit
48 200 03	10
48 200 04	10
48 200 05	10
48 200 06	20

Erlenmeyerkolben · Erlenmeyer flasks

Erlenmeyerkolben, Weithals, Polypropylen, transparent, mit Schraubkappe (NS-Stopfen lassen sich ebenfalls verwenden)

Erlenmeyer flasks, wide neck, polypropylene, translucent, with screw cap (NS stoppers can also be used)

Cat. No.	Inhalt Capacity	Gewinde Thread	für NS for NS	unit
58 110 00	50 : 10 ml	GL40	34/35	6
58 110 01	100 : 20 ml	GL40	34/35	6
58 110 02	250 : 50 ml	GL52	45/40	6
58 110 03	500 : 100 ml	GL52	45/40	6
58 110 04	1000 : 200 ml	GL52	45/40	4

VITLAB
Competence in Labware

