

Tracheal Cannulae

Tracheal Cannulae



Tracheal Cannula Only



Tracheal Cannula with Y-Adapter



Tracheal Cannula with Luer-Adapter

The tracheal cannula is used for general ventilation or respiratory mechanics studies. It is used in the case of a tracheotomy for acute studies. The outside diameter of the cannula must ensure a tight fit to the inside diameter of the animal's trachea. The cannulae are available in 3 versions:

- cannula only, for your own assembly
- with Y-adapter for connecting to a ventilator, the outside diameter of the Y-adapter matches to most of the ports diameter of the ventilators.
- with Luer-adapter for connecting to ventilation head of the multivent, see page F9.

Tracheotomy Cannulae (Stainless Steel)					
OD (mm)	Length [§] (mm)	Animal	Cannula only	with Y-adapter	with Luer-adapter
1.0 \$	20	mouse	BS4 73-0029	BS4 73-2731*	BS4 73-2835
1.2 \$	20	mouse	BS4 73-2825	BS4 73-2830*	BS4 73-2836
1.3 \$	20	mouse	BS4 73-0028	BS4 73-2730*	BS4 73-2837
1.5 \$	25	small rat	BS4 73-2826	BS4 73-2831 [#]	BS4 73-2729
1.8 \$	25	small rat	BS4 73-2827	BS4 73-2832 [#]	BS4 73-2728
2.0 \$	25	small rat	BS4 73-2828	BS4 73-2833 [#]	BS4 73-2727
2.3 \$	25	rat/guinea pig	BS4 73-2829	BS4 73-2834 [#]	BS4 73-2726
2.5 \$	25	rat/guinea pig	BS4 73-0033	BS4 73-2732 [#]	BS4 73-2725
3.0 \$	25	rat/guinea pig	BS4 73-0034	BS4 73-2733 [#]	BS4 73-2724
3.5 \$	30	rat/guinea pig	BS4 73-0035	BS4 73-2734 [#]	BS4 73-2723
4.0 \$	35	rabbit	BS4 73-0036	BS4 73-2735 [~]	-
5.0 \$	35	rabbit	BS4 73-0037	BS4 73-2736 [~]	-
6.0 \$	40	rabbit	BS4 73-3319	BS4 73-2891 [~]	-

* Y-adapter for Minivent OD 3.0 mm
[#] Medium Y-adapter OD 7.5 mm
[~] Large Y-adapter OD 10.0 mm
[§] Length listed is total cannula length. If cannula is used with Luer Adapter, the length is 4 mm less than listed. If the cannula is used with a Medium Y-adapter, the length is 3 mm less than listed. If the cannula is used with a Large Y-adapter, the length is 6 mm less than listed.

Intubation Cannulae



Intubation Cannula Only



Intubation Cannula with Y-Adapter




Intubation Cannula with Luer-Adapter

The intubation cannula is used in chronic studies. The outside diameter of the cannula must ensure a tight fit to the animal's trachea. These cannulae are also available in the 3 versions.

Intubation Cannulae (Material Stainless Steel)					
OD (mm)	Length [§] (mm)	Animal	Cannula only	with Y-adapter	with Luer-adapter
1.0 \$	28	mouse	BS4 73-0030	BS4 73-2737*	BS4 73-2848
1.1 \$	28	mouse	BS4 73-2838	BS4 73-2843*	BS4 73-2849
1.2 \$	30	mouse	BS4 73-2839	BS4 73-2844*	BS4 73-2850
1.5 \$	35	small rat	BS4 73-2840	BS4 73-2845 [#]	BS4 73-2851
1.8 \$	35	small rat	BS4 73-2841	BS4 73-2846 [#]	BS4 73-2852
2.0 \$	45	small rat	BS4 73-0038	BS4 73-2738 [#]	BS4 73-2853
2.3 \$	45	rat/guinea pig	BS4 73-2842	BS4 73-2847 [#]	BS4 73-2854
2.5 \$	50	rat/guinea pig	BS4 73-0039	BS4 73-2739 [#]	BS4 73-2855
3.0 \$	55	rat/guinea pig	BS4 73-0040	BS4 73-2740 [#]	BS4 73-2856
3.5 \$	60	rat/guinea pig	BS4 73-0041	BS4 73-2741 [~]	BS4 73-3112
4.5 \$	70	rabbit	BS4 73-0042	BS4 73-2742 [~]	-

* Y-adapter for Minivent OD 3.0 mm
[#] Medium Y-adapter OD 7.5 mm
[~] Large Y-adapter OD 10.0 mm
[§] Length listed is total cannula length. If cannula is used with Luer Adapter, the length is 4 mm less than listed. If the cannula is used with a Medium Y-adapter, the length is 3 mm less than listed. If the cannula is used with a Large Y-adapter, the length is 6 mm less than listed.



Catalog No. \$ **Product**
 BS4 73-2943 Adapter for Cannulae with Luer Connection and Pressure Measurement Port for KTR-4 Ventilator, pkg. of 5, see page F10

Replacement Parts

BS4 73-0027 Small Y Adapters, each with 2x silicone tubing 1.5 mm ID, 3 mm OD, 14 cm long, pkg. of 10
 BS4 73-0032 5 Medium Y adapters, each with 2x silicone tubing 6 mm ID, 9 mm OD, 80 cm long, pkg. of 5

Connect With Us



Previous



Next

Table of Contents

Product Index

Search

Home

Contact Us



Previous



Next



سوماتكو
 SOMATCO
 www.somatco.com