



W45110



W45151



W45025

Childbirth Simulator

This versatile simulator features the pelvic cavity with major anatomical landmarks and painted outline of bony pelvis. It allows intra uterine manipulation and various training possibilities for prenatal, intra partum, and postpartum skills as follows:

- Foetal palpation of fontanelle, backbone, knees, and elbows
- Vertex presentation
- Normal vaginal delivery
- Complete, frank and footling breech delivery
- C-section delivery
- Vertex and breech presentations in multiple delivery
- Demonstration of placenta previa: total, partial, and marginal
- Prolapse of umbilical cord
- Ritgen's manoeuvre
- 3 soft vulva inserts for episiotomy exercises
- Normal delivery of umbilical cord and placenta
- Suction of baby's nose and mouth

Supplied with 2 stomach covers (skin-tone /transparent), 1 male and 1 female foetus with placenta, 5 umbilical cords, 2 umbilical clamps, talcum powder and carrying bag.

53x33x43 cm; 8.0 kg

E

W45025



W45150

Options and Replacements for W45025

Articulating Foetus

Representation of a 42 cm foetus with articulating head and neck, shoulders, elbows, hips and knees for more realistic practice of difficult delivery exercises such as Leopold's or Ritgen's manoeuvres.

Features include:

- Realistic mouth and nostrils for simulated suctioning
- Palpable fontanelle, spine, shoulders, buttocks, elbows, and knees

Supplied with carrying bag.

81x20x20 cm; 2.3 kg

E

W45110

Palpation Module for Leopold's Manoeuvres

This module consists of a foetus, 2 cushions which may be inflated with air and an inflation device. Inflation of the lower cushion raises the foetus to desired position; inflation of the upper cushion creates a firm abdomen as in the ninth month of pregnancy. The foetus with palpable fontanelle, spine, shoulders, elbows, and knees may be placed in normal, breech, or transverse position and will be cradled between the 2 cushions.

33x30x18 cm; 1.8 kg

E

W45150

Labour Delivery Module

The 6 labour stations represent conditions of the cervix and vagina prior to labour, during labour, and at birth in a primigravid woman.

Stations illustrated are:

- Station -5: Prior to onset of labour
- Station -4: Cervix partially effaced
- Station -3: Cervix fully effaced
- Station 0: Foetal head at plane of ischial spine
- Station +2: Cervix nearing full dilation
- Station +5: Crowning of foetal head

Delivered with carrying case.

52x34.3x16.5 cm; 0.5 kg

E

W45151

Foetus for Vacuum Delivery (not shown)

Features smooth, soft skull with fontanelle for realistic vacuum delivery. With placenta and umbilical cord.

52x14x11.4 cm; 1.0 kg

E

W45026

Connect With Us

