



### **Gynaecology Simulator**

This simulator represents a female adult lower body with removable soft outer skin as well as anatomically accurate fimbriae and ovaries. Various pelvic examinations including recognition of pathologies and placing of contraceptive devices can be trained. Features include:

- Vaginal speculum and bi-manual pelvic examination
- · Uterine sounding
- · Insertion and removal of contraceptive devices such as IUD, diaphragm, cervical cap, contraceptive sponge or female condom

The simulator is supplied with normal and abnormal cervices and uteri which can be easily changed:

- 1 normal anteverted uterus with transparent top and round ligaments
- 1 normal uterus
- 6 uteri with external pathologies
- 2 normal cervices with patent as for IUD insertion/removal
- 6 cervices with pathologies

Supplied with talcum powder and carrying bag. Speculum is not included in scope of delivery.

51x46x25.5 cm; 8.5 kg

 $\square$  E

W45024

# **Options and Replacements for W45043**

7 Uteri - Normal and with External Pathologies  $\square$  E

W45044

7 Uteri - Normal and with Internal Pathologies For hysteroscopic viewing.

□ E

W45045

## **Postpartum Uterus**

48 hour postpartum uterus with duckbill cervix and fallopian tubes for IUD insertion, using long curved forceps. Speculum and bent forceps are not included in scope of delivery. 29.2x14x14 cm; 0.5 kg

W45101

### 2 Early Pregnant Uteri

6-8 week and 10-12 week pregnancy uteri, both with 3 interchangeable cervices.

W45102



#### **Gynaecology Simulator**

This simulator represents a female adult lower body with relevant internal anatomical landmarks and anatomically accurate fimbriae and ovaries. Various pelvic examinations including laparoscopy, mini laparotomy as well as recognition of pathologies and placing of contraceptive devices can be trained. Features include:

- Vaginal speculum and bi-manual pelvic examination
- Uterine sounding
- Insertion and removal of contraceptive devices such as IUD, 75mm diaphragm, or female condom
- Laparoscopic visualisation and occlusion of fallopian tubes
- Mini laparotomy

The simulator is supplied with normal and abnormal cervices as well as normal and pregnant uteri which can be easily changed:

- 1 anteverted and 1 retroverted parous uterus
- 2 pregnant uteri (ten week / twenty week)
- 5 normal cervices with patent as for IUD insertion / removal
- · 4 cervices with pathologies


Supplied with 10 fallopian tubes, talcum powder and carrying bag. The image shows an opened model to make the inner structures visible. The model is supplied with a closed abdomen. Hook is not included in scope of delivery.

51x46x25.5 cm; 8.5 kg

₽ E

W45043





Connect With Us











