



XC-401

XC - 401

MULTI FUNCTIONAL PATIENT CARE MANIKIN

The model is made of soft and half-hard PVC plastic, solid and durable, easy to use. All joints are movable, can be bent, all parts can be disassembled. It can demonstrate following items:

- | | |
|--|---------------------------------------|
| 1. Washing face and bathing the bed-ridden patient. | 12. Abdomen puncture |
| 2. Care and cleaning of the mouth | 13. Medulla puncture |
| 3. Care and cleaning of trachea dissection | 14. Lumbar puncture |
| 4. Inhalation of oxygen (nasal cork or nasal catheter) | 15. Deltoid injection |
| 5. Nasal lavage | 16. Deltoid hypodermic injection |
| 6. Stomach lavage | 17. Vein injection |
| 7. Heart infection | 18. Vein puncture |
| 8. Heart resuscitating on chest | 19. Vein transfusion |
| 9. Chest puncture. | 20. Vein blood transfusion |
| 10. Liver puncture | 21. Urethral catheterization of women |
| 11. Kidney puncture | 22. Muscular injection on buttocks |



XC-409

XC - 409A

NEW STYLE NEWBORN BABY MODEL (MALE)

XC - 409B

NEW STYLE NEWBORN BABY MODEL (FEMALE)

The size and configuration of this model is almost as same as the normal newborn baby. The head, upper limbs, lower limbs and all joints are flexible. It can demonstrate the characteristic of the newborn baby and all kind of operations in clinic education: includes different parturition method, newborn baby bath etc. Made of soft PVC plastic. Washable. Length: 52 CM, head diameter: 34CM.



XC - 402

XC - 402

MODEL FOR COURSE OF DELIVERY

The model basing on the mechanism of normal labor, shows the cardinal movements of engagement: descent-flexion-internal rotation of the head-extension-rotation of the head-external rotation of the head-internal rotation of the shoulders at accouchement. Suitable for teaching the mechanism of normal labor in medical colleges and nursing schools.



XC-410

XC - 410

NEWBORN INTUBATION BABY MODEL

This model is used to demonstrate the newborn baby intubation and train the students to learn intubation skills. Served as an ideal model for medical college, nursing schools and hospital etc. Made of elastic PVC plastic.

XC - 414

DEVELOPMENT PROCESS FOR FETUS

1. Normal uterus, shows the inner genitals of female.
2. Shows the embryo and uterus in the first month of gestation, placenta is shown also.
3. Shows the fetus & uterus in the second month of gestation, the embryo is just looks like human.
4. Shows the fetus and uterus in the third month of gestation, the height of the fetus is about 9 cm.
5. Shows the fetus and uterus in the fourth month of gestation, the height of the fetus is about 16 cm.
6. Shows the fetus and uterus in the fifth month of gestation, the height of the fetus is about 25 cm.
7. Shows the fetus and uterus in the sixth month of gestation, the height of the fetus is about 30 cm.
8. Shows the fetus and uterus in the seventh month of gestation, the height of the fetus is about 35 cm.
9. Shows the fetus and uterus in the eighth month of gestation, the height of the fetus is about 40 cm.
10. Shows the fetus and uterus in the ninth month of gestation, the height of the fetus is about 45 cm.

Dimension: Natural size, on the base.



XC-414

Connect With Us

